# EMP 4 ECO+

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## Content

|                                          |    |
|------------------------------------------|----|
| General Information .....                | 38 |
| Description of the EMP4 ECO +.....       | 39 |
| Description of the Symbols .....         | 39 |
| Technical Specifications.....            | 40 |
| Operating the EMP4 ECO + .....           | 41 |
| Setting of Intensity Levels.....         | 45 |
| Principle of Dynamic Stimulation .....   | 46 |
| Short Instruction.....                   | 47 |
| Description of the Pre-Set Programs..... | 48 |
| Technical Information .....              | 58 |
| Accessories .....                        | 60 |
| Delivery Content .....                   | 62 |
| Electrodes Position.....                 | 63 |
| Program Overview .....                   | 70 |

## General Information

### Intended Use

The EMP4 ECO + device has been designed for the generation of transcutaneous stimulation of muscles and nerves in humans and may not be used for any other purposes.

### Safety Instructions

Please read the operating instructions carefully before using the EMP4 ECO +.

#### **To ensure safe use of the EMP4 ECO +:**

Only use the EMP4 ECO + with original accessories.

Keep the EMP4 ECO + away from water or other liquids.

Do not drop the EMP4 ECO +, do not use it inappropriately or expose it to extreme temperatures or high levels of humidity (not less than 10°C or more than 40 °C or a relative humidity of more than 90 %).

Do not use the EMP4 ECO + if it is not working properly or if it has been damaged in any way.

Be careful when using the EMP4 ECO + near or on children.

Always store the EMP4 ECO + in its case to protect it from damage and dust.


The size of the electrodes are to be not less than 2 cm<sup>2</sup>.

**Caution !** Patients with medical implants (e.g. cardiac pacemakers) should consult their doctor before using the EMP4 ECO +. The EMP4 ECO + may only be connected to one person at any one time.


Long treatment sessions are subject to limited periods of time. Do not extend the treatment time of the fixed programs. After a treatment session, wait for several hours before using the EMP4 ECO + again.

**Warning !** The simultaneous use of EMP4 ECO + and a high frequency surgery unit can cause burns underneath the electrodes. The application of the unit close to a shortwave or microwave device can cause variations in the output values of the EMP4 ECO +. In addition, portable or mobile high frequency communication devices may affect the operation of electrical medical devices.

## Description of the EMP4 ECO +

1. Display
2. Program selection key
3. To enter the user mode and to pause the stimulation
4. Keys to adjust the intensities of all channels
  - ⤴ Increasing intensity
  - ⤵ Decreasing intensity
5. ON/OFF key
6. Output sockets for the cables
7. Output socket for the battery charger
8. ON/OFF commutator

## Description of the Symbols

Attention: Read accompanying documents, especially user manual!

By labelling with CE mark, the manufacturer asserts that the product fulfills all relevant requirements of the respective EC Directives. A conformation assessment process has been successfully completed. The identification number in respect of the conformation assessment process is given in accordance with the CE labelling of the relevant notified body.

Type BF applied part

Serial Number

Catalogue number

Date of manufacture

Manufacturer

Environment

Do not dispose of the appliance with the normal household waste at the end of its life, but hand it in at an official collection point for recycling. By doing this you will help to protect the environment.

## Technical Specifications


Nerve and muscle stimulator with 2x2 channels (4x100 mA max.) electrically insulated channels, constant current characteristic, 25 integrated pre-set programs and 20 user programs.

|                     |                                                  |
|---------------------|--------------------------------------------------|
| Output current      | 100 mA (with an actual resistor of 1k $\Omega$ ) |
| Frequency range     | 1-120 Hz                                         |
| Impulse width       | 50-500 $\mu$ s                                   |
| Current consumption | 15 mA                                            |
| Power supply        | Rechargeable battery Ni-MH 4.8 V                 |
| Dimensions          | 120 x 75 x 35 mm                                 |
| Weight              | approx. 250 g                                    |

**Impulse form:** Biphasic symmetric impulses with pause

With 1 k $\Omega$  load real

With ANSI standard load

## Operating the EMP4 ECO +


### Connecting the Cables and the Electrodes

Connect the electrodes with the cable(s) (always 2 electrodes per cable). Connect the cables to the EMP4 ECO + on sockets no. 6. Place the electrodes on the desired area (see electrode placement)

### Switching ON the EMP4 ECO +


Switch the commutator (no. 8) to the ON position. Press the  key to start the device. The program shown on the display is the last program used. The device automatically switches off after 2 minutes of disuse. If the display shows a key symbol on the upper left side, the EMP4 ECO + has previously locked a specific program (see chapter locking the key pad).

### Selecting a Program

Ensure that the EMP4 ECO + has not been locked on one program. Select a program by pressing the **P** key. To scroll through the programs, press the **P** key as often as necessary until you have selected the program of your choice. To return to the previous program, press one of the  keys.

### Starting the Stimulation

First check if you have selected the right program and if cables and electrodes are well connected and well positioned on the skin. Start the stimulation by increasing the intensity with the  key of the channel(s) connected to the electrodes.

### Setting the Intensity

Adjust the intensity of the channels  (black) and  (red) by using the left  and  keys.

Adjust the intensity of the channels  (grey) and  (yellow) by using the right  and  keys.

Increase the intensity with the  keys.

Decrease the intensity with the  keys.

By keeping the intensity keys pressed, the intensity will increase/decrease more rapidly. The intensity of each channel is displayed on the screen. The intensity can be set between 0 and 100 mA.

**Caution!** The stimulator automatically locks the intensity to the set level after 5 seconds. To adjust the intensity you firstly need to decrease the intensity. The intensity can then be increased again.

If electrodes are not or not correctly connected to the stimulator, you will not be able to set the intensity to more than 10 mA. The same result will occur if electrodes are old or the cables broken.

### Stopping the Stimulation

You can stop the stimulation at any time by pressing the **P** key or the  key for 2 seconds. When the therapy time is completed, the stimulation stops automatically.

### Pausing the Stimulation

You can temporarily stop the stimulation by pressing the **E** key during the stimulation. The word "PAUSE" then flashes on the display. You can then adjust the position of the electrodes. Press the E key again to re-start the stimulation. The intensity levels slowly increase to the level set before the pause.

## Creating User Programs (U1 to U20 only)

Choose one program between U1 and U20. Press the E key to access the Programming mode. The parameters to be changed are flashing. You can modify them by pressing the  and  keys. Switch to the next parameter by pressing the **E** key. The frequency, the pulse width and the time can be modified in the free programs.


The parameters of the user programs can be modified as follows:

|                                                                                                                                                                                                           |                          |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|
| U1, U3, U4, U9, U20                                                                                                                                                                                       |                          |
| Frequency                                                                                                                                                                                                 | 20-120 Hz                |
| Pulse width                                                                                                                                                                                               | 70-500 $\mu$ s           |
| Time                                                                                                                                                                                                      | 1-99 min                 |
| U2, U5 (TENS)                                                                                                                                                                                             |                          |
| Frequency channel  (black) and  (red)   | 80, 90, 100, 110, 120 Hz |
| Frequency channel  (grey) and  (yellow) | 2, 5, 10 Hz              |
| Pulse width                                                                                                                                                                                               | 70-500 $\mu$ s           |
| Time                                                                                                                                                                                                      | 1-99 min                 |
| U6                                                                                                                                                                                                        |                          |
| Time                                                                                                                                                                                                      | 1-99 min                 |
| U7, U8                                                                                                                                                                                                    |                          |
| Maximum Frequency                                                                                                                                                                                         | 10-120 Hz                |
| Minimum Frequency                                                                                                                                                                                         | 2 Hz                     |
| Pulse width                                                                                                                                                                                               | 70-180 $\mu$ s           |
| Time                                                                                                                                                                                                      | 1-99 min                 |
| U10, U11, U17                                                                                                                                                                                             |                          |
| Frequency                                                                                                                                                                                                 | 1-120 Hz                 |
| Pulse width                                                                                                                                                                                               | 70-500 $\mu$ s           |
| Time                                                                                                                                                                                                      | 1-99 min                 |
| U12, U13, U14, U15, U16, U18, U19                                                                                                                                                                         |                          |
| Frequency                                                                                                                                                                                                 | 20-120 Hz                |
| Pulse width                                                                                                                                                                                               | 70-500 $\mu$ s           |
| Rising ramp                                                                                                                                                                                               | 1-4 s                    |
| Working time                                                                                                                                                                                              | 1-25 s                   |
| Pause                                                                                                                                                                                                     | 1-25 s                   |
| Time                                                                                                                                                                                                      | 1-99 min                 |

## Switching-OFF the Stimulator

Press the  key for one second to stop the stimulation and press  key a second time to turn OFF the EMP 4 ECO +. If the battery level is too low or if the stimulator has not been used for the previous 2 minutes, it turns off automatically.

## Charging the Battery

Output socket  
for the battery  
charger

The battery level is shown on the display with a battery symbol with 4 lines. If the battery level is too low, the stimulator switches off automatically and cannot be turned on again. The battery then needs to be re-charged. Only use the delivered charger (Item no. 101062) to charge the battery of the EMP4 ECO +.

Switch the commutator to the OFF position.

Connect the charger with the EMP4 ECO +.

Connect the charger to the mains circuit. The diode on the charger illuminates in red.


Charge the EMP4 ECO + until the diode on the charger illuminates in green.

The battery is fully charged if the light on the charger appears in green.

Disconnect the charger from the main circuit and from the stimulator. Place the charger back in the storage case.

Important! Do not forget to switch the commutator to ON again after charging.

Charger

## Setting of Intensity Levels

Do not try to achieve a higher and higher level of intensity. Set the intensity in such a way that you get a comfortable feeling during the stimulation. Increase the intensity carefully to a maximum tolerable and then reduce it slightly step by step to a comfortable level. Note that the intensity settings change depending on the stimulation area and the time. This is normal and can be explained as follows:

**Resistance of the skin:** Dry skin has a lower conductivity level than damp skin (due to sweating). The skin does not show the same resistance to linear current on all parts of the body. For example, the resistance level of horny skin can be twice that of the skin in the hollow of the knee. The innervation differs greatly according to the area of the body. This results in a difference of sensitivity.

**Volume of the muscle mass:** If you stimulate muscles with a high volume you will automatically need to set a higher intensity than for the stimulation of small muscles.

**Muscle fatigue:** A tired muscle will only tolerate a lower intensity.

**Age of electrodes:** Self-adhesive electrodes do not last forever. You can use them up to for up to 60 stimulation sessions. This depends on the quality of the electrodes but also on the skin of the user. Old electrodes have a higher resistance and their conductivity levels decline. It is important to replace electrodes regularly.

**Adaption to the current:** The nerves adapt rapidly to the current. It is normal to feel "less current" after a few minutes of stimulation with specific parameters. In such a case feel free to increase the intensity in order to get the best results.


**Frequency and pulse width:** Frequencies and pulse widths differ depending on the programs. Therefore the intensity cannot be set at the same lever for all the programs even if the electrodes are placed in the same position. The higher the frequency is set, the lower you will be able to set the intensity.

The higher the pulse width is set, the lower you will be able to set the intensity.

## Principle of Dynamic Stimulation

The dynamic stimulation principle is based on the fact that the stimulation is not applied simultaneously on all channels. Dynamic stimulation generates a wave movement of the stimulation impulse which moves to and fro between the electrode pairs and therefore attains enhanced reproduction of the massaging effect. Dynamic stimulation creates for the user a more pleasurable experience than with classic neuro-muscular stimulation. The consecutive four channel distribution of the electrical current produces good results in respect of both pain alleviation and muscular relaxation. The application of dynamic stimulation for lymph drainage accentuates the effect of the wave type muscular contraction.

### Electrode Positions for dynamic Stimulation

Dynamic stimulation on legs: Place one electrode of each channel on each leg (size 5x9 cm). All anodes should be on one leg and all cathodes on the other leg.

Dynamic stimulation on the back: Place one electrode of each channel on each side of the spinal column (size 5x9 cm). Place all anodes on one side and all cathodes on the other side.

### Important note

If you use electrodes with two connectors, always connect them to the same pole colour of the cables (either both blue or both red).

## Short Instruction

### EMP 4 ECO+ SHORT INSTRUCTION

**Important !** Connect the electrodes with the cable(s) (always 2 electrodes per cable). Connect the cables to the device. Place the electrodes on the desired area.

#### 1 Switch ON/OFF the unit

#### 2 Program choice

Press the **P** key to select a program.

#### 3 Intensity setting

Intensity channel 1

increase the intensity  
decrease the intensity


Intensity channel 2

increase the intensity  
decrease the intensity

**Caution!** The intensity is automatically interlocked after 5 s. First decrease the intensity to increase again the intensity.

#### 4 Stop stimulation


Press the **P** key for 2 seconds to stop the stimulation.

# Description of the Pre-Set Programs


## Program P1 + U1:

### Low Back Pain

|                                           |                                                                                                                                                                                                   |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Type of stimulation                       | Gate control 80 Hz - dynamic stimulation (slow)                                                                                                                                                   |
| Parameter                                 | Frequency: 80 Hz                                                                                                                                                                                  |
|                                           | Pulse width: 150 µs                                                                                                                                                                               |
|                                           | Time: 20 min                                                                                                                                                                                      |
|                                           | Rising ramp: 1 s                                                                                                                                                                                  |
|                                           | Falling ramp: 1 s                                                                                                                                                                                 |
| Symbols on the display during stimulation | Rising phase <br>Falling phase  |
| Description                               | The 4 channels are working with the same frequency and pulse width.<br>                                          |

## Program P2 + U2:

### Lumbosciatica

|                                          |                                                                                                                                                                                                      |
|------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Stimulationsart                          | Gate Control 80 Hz + 2 Hz Endorphinausschüttung                                                                                                                                                      |
| Parameter                                | Frequency Channel  and  : 80 Hz  |
|                                          | Frequency Channel  and  : 2 Hz |
|                                          | Pulse width: 200 µs                                                                                                                                                                                  |
|                                          | Time: 20 min                                                                                                                                                                                         |
| Symbol on the display during stimulation |                                                                                                                   |
| Description                              |                                                                                                                   |

**Program P3 + U3:****Cervicobrachial Syndrome**

|                     |                                                 |
|---------------------|-------------------------------------------------|
| Type of stimulation | Gate control 80 Hz - dynamic stimulation (fast) |
|---------------------|-------------------------------------------------|

|           |                  |
|-----------|------------------|
| Parameter | Frequency: 80 Hz |
|-----------|------------------|

|  |                          |
|--|--------------------------|
|  | Pulse width: 150 $\mu$ s |
|--|--------------------------|

|  |              |
|--|--------------|
|  | Time: 20 min |
|--|--------------|

|  |                    |
|--|--------------------|
|  | Rising ramp: 0.2 s |
|--|--------------------|

|  |                     |
|--|---------------------|
|  | Falling ramp: 0.2 s |
|--|---------------------|


|                                           |  |
|-------------------------------------------|--|
| Symbols on the display during stimulation |  |
|-------------------------------------------|--|

|  |                                                                                                |
|--|------------------------------------------------------------------------------------------------|
|  | Rising phase  |
|--|------------------------------------------------------------------------------------------------|

|  |                                                                                                 |
|--|-------------------------------------------------------------------------------------------------|
|  | Falling phase  |
|--|-------------------------------------------------------------------------------------------------|

|             |  |
|-------------|--|
| Description |  |
|-------------|--|

|  |                                                                     |
|--|---------------------------------------------------------------------|
|  | The 4 channels are working with the same frequency and pulse width. |
|--|---------------------------------------------------------------------|

**Program 4 + U4:****Chronic Pain**

|                     |                              |
|---------------------|------------------------------|
| Type of stimulation | Burst 2 Hz endorphin release |
|---------------------|------------------------------|

|           |                   |
|-----------|-------------------|
| Parameter | Frequency: 100 Hz |
|-----------|-------------------|

|  |                          |
|--|--------------------------|
|  | Pulse width: 150 $\mu$ s |
|--|--------------------------|

|  |              |
|--|--------------|
|  | Time: 20 min |
|--|--------------|

|  |                       |
|--|-----------------------|
|  | Working phase: 0.25 s |
|--|-----------------------|


|  |               |
|--|---------------|
|  | Pause: 0.25 s |
|--|---------------|

|                                            |  |
|--------------------------------------------|--|
| Symbols on the display during stimulation: |  |
|--------------------------------------------|--|

|  |                                                                                                   |
|--|---------------------------------------------------------------------------------------------------|
|  | Working phase  |
|--|---------------------------------------------------------------------------------------------------|


|  |                                                                                           |
|--|-------------------------------------------------------------------------------------------|
|  | Pause  |
|--|-------------------------------------------------------------------------------------------|

|             |  |
|-------------|--|
| Description |  |
|-------------|--|


Program 5 + U5:

Gate Control + Endorphin Release

|                                          |                                                                                                                                                                                                                                                                                      |
|------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Type of stimulation                      | Gate control on channel 1 + endorphin release on channel 2                                                                                                                                                                                                                           |
| Parameter                                | Frequency channel ● (black) and ● (red)(red): 100 Hz<br>Frequency channel ● (grey) and ● (yellow): 2 Hz<br>Pulse width: 200 µs<br>Time: 20 min                                                                                                                                       |
| Symbol on the display during stimulation | ■                                                                                                                                                                                                                                                                                    |
| Description                              | Channel ● (black) and ● (red) are working with a 100 Hz frequency (80-120 Hz in the user mode). Channel ● (grey) and ● (yellow) are working with a 2 Hz frequency (0.5-10 Hz in the user mode).<br> |

Program 6 + U6:


Gonarthrosis / Coxarthrosis

|                        |                                                                                                                          |
|------------------------|--------------------------------------------------------------------------------------------------------------------------|
| Type of stimulation    | Han                                                                                                                      |
| Parameter phase 1      | Frequency: 100 Hz<br>Pulse width: 100 µs<br>Time: 3 s                                                                    |
| Parameter phase 2      | Frequency: 2 Hz<br>Pulse width: 200 µs<br>Time: 3 s<br>Time: 20 min Phase 1 and phase 2 are alternating every 3 seconds. |
| Symbols on the display | Phase 1 ■<br>Phase 2 ■                                                                                                   |
| Description            |                                       |

## Program 7 + U7:

## Epicondylitis


|                                           |                                                                                                                                                                                                   |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Type of stimulation                       | Modulation                                                                                                                                                                                        |
| Parameter                                 | Minimum frequency: 2 Hz (Pulse width: 200 $\mu$ s)                                                                                                                                                |
|                                           | Maximum frequency: 80 Hz (Pulse width: 100 $\mu$ s)                                                                                                                                               |
|                                           | Duration of the modulation: 7.5 s                                                                                                                                                                 |
|                                           | Time: 20 min                                                                                                                                                                                      |
|                                           | Rising ramp: 0,2 s                                                                                                                                                                                |
|                                           | Falling ramp: 0,2 s                                                                                                                                                                               |
| Symbols on the display during stimulation | Rising phase <br>Falling phase  |
| Description                               |                                                                                                                                                                                                   |

## Program 8 + U8:


## Algodystrophy

|                                           |                                                                                                                                                                                                       |
|-------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Type of stimulation                       | Modulation                                                                                                                                                                                            |
| Parameter                                 | Minimum frequency: 2 Hz (Pulse width: 148 $\mu$ s)                                                                                                                                                    |
|                                           | Maximum frequency: 80 Hz (Pulse width: 70 $\mu$ s)                                                                                                                                                    |
|                                           | Duration of the modulation: 7.5 s                                                                                                                                                                     |
|                                           | Time: 20 min                                                                                                                                                                                          |
|                                           | Rising ramp: 0.2 s                                                                                                                                                                                    |
|                                           | Falling ramp: 0.2 s                                                                                                                                                                                   |
| Symbols on the display during stimulation | Rising phase <br>Falling phase  |
| Description                               |                                                                                                                                                                                                       |


Program 9 + U9:

Anti-Adaptation – Dynamic Stimulation

|                                           |                                                                                                                                                                                                   |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Type of stimulation                       | Dynamic stimulation (fast)                                                                                                                                                                        |
| Parameter                                 | Frequency: 80 Hz                                                                                                                                                                                  |
|                                           | Pulse width: 150 µs                                                                                                                                                                               |
|                                           | Time: 20 min                                                                                                                                                                                      |
|                                           | Rising ramp: 0.5 s                                                                                                                                                                                |
|                                           | Falling ramp: 0.5 s                                                                                                                                                                               |
| Symbols on the display during stimulation | Rising phase <br>Falling phase  |
| Description                               | The 4 channels are working with the same frequency and pulse width.<br>                                          |

Program 10 + U10:

Classic Gate Control 120 Hz

|                                          |                                                                                     |
|------------------------------------------|-------------------------------------------------------------------------------------|
| Type of stimulation                      | TENS                                                                                |
| Parameter                                | Frequency: 120 Hz                                                                   |
|                                          | Pulse width: 150 µs                                                                 |
|                                          | Time: 20 min                                                                        |
|                                          | Working phase: 0.25 s                                                               |
|                                          | Pause: 0.25 s                                                                       |
| Symbol on the display during stimulation |  |
| Description                              |  |

**Program 11 + U11:****Treatment of Contractures**


|                                          |                          |
|------------------------------------------|--------------------------|
| Type of stimulation                      | TENS                     |
| Parameter                                | Frequency: 1 Hz          |
|                                          | Pulse width: 150 $\mu$ s |
|                                          | Time: 20 min             |
|                                          | Working phase: 0.25 s    |
|                                          | Pause: 0.25 s            |
| Symbol on the display during stimulation | ■                        |
| Description                              | See picture program P10  |

**Muscle Stimulation Programs**

|                        |                 |
|------------------------|-----------------|
| Symbols on the display | Rising phase ▲  |
|                        | Working phase ■ |
|                        | Pause —         |

**Program 12 + U12:****Atrophy Upper Extremities**

|                     |                          |
|---------------------|--------------------------|
| Type of stimulation | Muscle stimulation       |
| Parameter           | Frequency: 35 Hz         |
|                     | Pulse width: 150 $\mu$ s |
|                     | Rising ramp: 2 s         |
|                     | Working phase: 3 s       |
|                     | Pause: 9 s               |
|                     | Time: 20 min             |

**Program 13 + U13:****Atrophy Lower Extremities**

|                     |                          |
|---------------------|--------------------------|
| Type of stimulation | Muscle stimulation       |
| Parameter           | Frequency: 35 Hz         |
|                     | Pulse width: 300 $\mu$ s |
|                     | Rising ramp: 2 s         |
|                     | Working phase: 3 s       |
|                     | Pause: 9 s               |
|                     | Time: 20 min             |
| Description         | See picture program P12  |

**Program 14 + U14:****Muscle Strengthening Upper Extremities**

|                     |                         |
|---------------------|-------------------------|
| Type of stimulation | Muscle stimulation      |
| Parameter           | Frequency: 65 Hz        |
|                     | Pulse width: 250 µs     |
|                     | Rising ramp: 2 s        |
|                     | Working phase: 4 s      |
|                     | Pause: 8 s              |
|                     | Time: 20 min            |
| Description         | See picture program P12 |

**Programm 15 + U15:****Muscle Strengthening Lower Extremities**

|                     |                         |
|---------------------|-------------------------|
| Type of stimulation | Muscle stimulation      |
| Parameter           | Frequency: 65 Hz        |
|                     | Pulse width: 300 µs     |
|                     | Rising ramp: 2 s        |
|                     | Working phase: 4 s      |
|                     | Pause: 8 s              |
|                     | Time: 20 min            |
| Description         | See picture program P12 |

**Programm 16 + U16:****Venous Reflux**

|                     |                         |
|---------------------|-------------------------|
| Type of stimulation | Muscle stimulation      |
| Parameter           | Frequency: 35 Hz        |
|                     | Pulse width: 250 µs     |
|                     | Rising ramp: 3 s        |
|                     | Working phase: 5 s      |
|                     | Pause: 10 s             |
|                     | Time: 20 min            |
| Description         | See picture program P12 |

**Program 17 + U17:****Urge Incontinence**

|                     |                         |
|---------------------|-------------------------|
| Type of stimulation | TENS                    |
| Parameter           | Frequency: 10 Hz        |
|                     | Pulse width: 180 µs     |
|                     | Time: 15 min            |
| Description         | See picture program P10 |

**Program 18 + U18:****Mixed Incontinence**


|                     |                          |
|---------------------|--------------------------|
| Type of stimulation | Muscle stimulation       |
| Parameter           | Frequency: 20 Hz         |
|                     | Pulse width: 180 $\mu$ s |
|                     | Rising ramp: 2 s         |
|                     | Working phase: 4 s       |
|                     | Pause: 4 s               |
|                     | Time: 15 min             |
| Description         | See picture program P12  |

**Program 19 + U19:****Stress Incontinence**

|                     |                          |
|---------------------|--------------------------|
| Type of stimulation | Muscle stimulation       |
| Parameter           | Frequency: 50 Hz         |
|                     | Pulse width: 180 $\mu$ s |
|                     | Rising ramp: 2 s         |
|                     | Working phase: 3 s       |
|                     | Pause: 6 s               |
|                     | Time: 15 min             |
| Description         | See picture program P12  |

**Program 20 + U20:****Relaxation – Dynamic Stimulation**

|                                           |                                                                                                                                                                                                       |
|-------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Type of stimulation                       | Dynamic stimulation (slow)                                                                                                                                                                            |
| Parameter                                 | Frequency: 120 Hz                                                                                                                                                                                     |
|                                           | Pulse width: 150 $\mu$ s                                                                                                                                                                              |
|                                           | Time: 20 min                                                                                                                                                                                          |
|                                           | Rising ramp: 2 s                                                                                                                                                                                      |
|                                           | Falling ramp: 2 s                                                                                                                                                                                     |
| Symbols on the display during stimulation | Rising phase <br>Falling phase  |
| Description                               | The 4 channels are working with the same frequency and pulse width.                                                                                                                                   |


Program 21 + U21:

Combination: Active wrist joint motion + pain management

|                     |                                                                                                                        |
|---------------------|------------------------------------------------------------------------------------------------------------------------|
| Type of stimulation | Muscle stimulation and gate control                                                                                    |
| Parameter           | Frequency channel <span style="color: black;">●</span> (black) and <span style="color: red;">●</span> (red): 45 Hz     |
|                     | Frequency channel <span style="color: grey;">●</span> (grey) and <span style="color: yellow;">●</span> (yellow): 90 Hz |
|                     | Pulse width: 200 µs                                                                                                    |
|                     | Rising ramp: 3 s                                                                                                       |
|                     | Working phase: 1 s                                                                                                     |
|                     | Falling ramp: 1 s                                                                                                      |
|                     | Pause: 6 s                                                                                                             |
|                     | Time: 20 min                                                                                                           |

Description

Program 22 + U22:

Combination: Active elbow joint motion + pain management

|                     |                                                                                                                        |
|---------------------|------------------------------------------------------------------------------------------------------------------------|
| Type of stimulation | Muscle stimulation and gate control                                                                                    |
| Parameter           | Frequency channel <span style="color: black;">●</span> (black) and <span style="color: red;">●</span> (red): 45 Hz     |
|                     | Frequency channel <span style="color: grey;">●</span> (grey) and <span style="color: yellow;">●</span> (yellow): 90 Hz |
|                     | Pulse width: 250 µs                                                                                                    |
|                     | Rising ramp: 5 s                                                                                                       |
|                     | Working phase: 1 s                                                                                                     |
|                     | Falling ramp: 5 s                                                                                                      |
|                     | Pause: 10 s                                                                                                            |
|                     | Time: 20 min                                                                                                           |

Description

See picture program P21

**Program 23 + U2 :****Combination: Active knee joint motion + pain management**

|                     |                                                  |
|---------------------|--------------------------------------------------|
| Type of stimulation | Muscle stimulation and gate control              |
| Parameter           | Frequency channel ● (black) and ● (red): 45 Hz   |
|                     | Frequency channel ● (grey) and ● (yellow): 90 Hz |
|                     | Pulse width: 300 µs                              |
|                     | Rising ramp: 5 s                                 |
|                     | Working phase: 1 s                               |
|                     | Falling ramp: 5 s                                |
|                     | Pause: 10 s                                      |
|                     | Time: 20 min                                     |
| Description         | See picture program P21                          |

**Program 24 + U24:****Combination: Active ankle joint motion + pain management**

|                     |                                                  |
|---------------------|--------------------------------------------------|
| Type of stimulation | Muscle stimulation and gate control              |
| Parameter           | Frequency channel ● (black) and ● (red): 45 Hz   |
|                     | Frequency channel ● (grey) and ● (yellow): 90 Hz |
|                     | Pulse width: 300 µs                              |
|                     | Rising ramp: 3 s                                 |
|                     | Working phase: 1 s                               |
|                     | Falling ramp: 1 s                                |
|                     | Pause: 6 s                                       |
|                     | Time: 20 min                                     |
| Description         | See picture program P21                          |

**Program 25 + U25:****Combination: Hemiparesis (shoulder) + pain management**

|                     |                                                  |
|---------------------|--------------------------------------------------|
| Type of stimulation | Muscle stimulation and gate control              |
| Parameter           | Frequency channel ● (black) and ● (red): 45 Hz   |
|                     | Frequency channel ● (grey) and ● (yellow): 90 Hz |
|                     | Pulse width: 250 µs                              |
|                     | Rising ramp: 2 s                                 |
|                     | Working phase: 6 s                               |
|                     | Falling ramp: 2 s                                |
|                     | Pause: 10 s                                      |
|                     | Time: 20 min                                     |
| Description         | See picture program P21                          |

## Technical Information

### Adjustments, Modifications and Repairs

The manufacturer is only responsible for the safety and performance of the EMP4 ECO + when readjustments, alterations and repairs are carried out by authorized persons and when the EMP4 ECO + is used in accordance with the operating instructions.

### Warranty

Legal right of warranty is applied according to German Civil Code.

### Guarantee

The manufacturer issues a guarantee of 12 months from the date of purchase.

The guarantee does not apply in the following cases:

- Damage due to improper handling
- Defects the customer is aware of on the date of purchase
- Damage caused by the customer
- For wearing parts and consumable supplies like, for instance, battery, cables and electrodes.

### Maintenance and Cleaning

No special cleaning or care agents are required for the EMP4 ECO +. Clean the stimulator with a soft, lint-free cloth. Please ensure that no moisture permeates the stimulator. If moisture does permeate the stimulator, a technical check must be carried out before re-use.

### Battery Retraction and Disposal

In connection with the sales of products which contain batteries or rechargeable batteries, we are obliged by law (battery law § 18, section 1) to advise you of the following:

The dustbin symbol () marks batteries with hazardous substances as well as the circumstance that batteries must not be disposed of in the household waste but professionally. The chemical term of the hazardous substance is below the dustbin symbol.

You are obliged by law to return used batteries.

You may render used batteries at municipal collecting points or in

local trade. We are also obliged to take back used batteries, whereat our obligation restricts to those which are or were part of our product range.

Therefore, you may send back used batteries to us adequately stamped by mail or render these at our warehouse free of charge at the following address:

schwa-medico GmbH, Wetzlarer Str. 41-43, 35630 Ehringshausen, Germany

Please refer to the following figures to see the respective symbols of batteries with hazardous substances:

Battery contains more than 0.002 % cadmium

Battery contains more than 0.0005 % quicksilver

Battery contains more than 0.004 % plumb

## Classification

The EMP4 ECO + has been designed as a class IIa product in accordance with appendix IX of the Medical Device Directive 93/42/EEC.

## Technical Check

We recommend a technical check on the EMP4 ECO + every 24 months. This includes:

1. Confirming that user instructions have been included in the accompanying documentation.
2. Checking the equipment for completeness.
3. Visual check:
  - for mechanical damage
  - for damage to all cables and plugs / sockets
4. Functional Safety
  - Checking the output signals with a load resistance of 1 kΩ real current and voltage
  - Checking the frequency
  - Checking the pulse width


These technical checks may only be performed by suitably qualified persons. The results must be documented with the date and name of the person carrying out the check.

Combination with  
Original Accessories

The EMP4 ECO + is a CE marked medical device and may therefore only be used together with the accessories listed in the chapter Accessories of this manual. The manufacturer cannot guarantee the security of the stimulator if not used with original accessories.

Accessories

STIMEX  
Self-Adhesive Electrodes

|               |                                 |
|---------------|---------------------------------|
| Compound:     | conductive and bonding material |
| Manufacturer: | Pierenkemper GmbH, Germany      |

**Application:** Place self-adhesive electrodes directly on the appropriate skin treatment area, no special preparation of electrodes is necessary. See „Electrode Maintenance“ below to ensure the electrode’s effectiveness. Do not apply on broken skin !



**Electrode maintenance:** For hygienic reasons, self-adhering electrodes are for single patient use only. After every treatment session return electrodes to their original foil and store them in their plastic bag. To increase longevity, store the STIMEX electrodes in a cool area (e.g. refrigerator).

| Item no. | Size                  | Quantity |
|----------|-----------------------|----------|
| 281000   | STIMEX, round 32 mm Ø | 4        |
| 282000   | STIMEX, round 50 mm Ø | 4        |
| 283400   | STIMEX, 50 x 50 mm    | 4        |
| 283600   | STIMEX, 50 x 90 mm    | 2        |
| 283000   | STIMEX, 50 x 130 mm   | 2        |
| 283100   | STIMEX, 80 x 130 mm   | 2        |


Probes and Gel

Vaginal and anal probe are only to be used with the incontinence programs of this unit. The probe may be used by one person only and in the interests of hygiene should not be shared with other persons. The probe is reusable if it has been cleaned and stored correctly after each treatment session. The probe should only be used on intact mucous membrane. Do not use the probe if it is damaged or if it does

not work correctly. Do not use the probe if the packaging is damaged. Clean the probe after each electrotherapy session with water and soap or clean it with 70% alcohol.

| Item no. | Description                                                                                             | Quantity |
|----------|---------------------------------------------------------------------------------------------------------|----------|
| 103025   | Vaginal probe V2B<br>  | 1        |
| 103045   | Anal probe 12 C<br>    | 1        |
| 103044   | Analprobe TL1R<br>     | 1        |
| 101509   | Sonogel for probes<br> | 1        |

## Cables and Charger

| Item no. | Description                                                                                                   | Quantity |
|----------|---------------------------------------------------------------------------------------------------------------|----------|
| 106351   | Cable type 5.15 (black)                                                                                       | 1        |
| 106352   | Cable type 5.16 (grey)                                                                                        | 1        |
| 106353   | Cable type 5.17 (red)                                                                                         | 1        |
| 106354   | Cable type 5.18 (yellow)                                                                                      | 1        |
|          |                            |          |
| 101062   | Charger for EMP4 ECO +<br> | 1        |

STIMEX-TEX  
Stimulation Sleeves

Instruction: Moisten the inner side of STIMEX-TEX before every utilisation. Connect the STIMEX-TEX with the cable of the AmpliMove medical. Slip the STIMEX-TEX over the arm, elbow, foot, leg or knee. Do not apply STIMEX-TEX on wounds. Air-dry the STIMEX-TEX after every utilisation (not above a radiator). STIMEX-TEX may be washed by hand with cold water.

| Item no. | Description                                                 | Quantity |
|----------|-------------------------------------------------------------|----------|
| 107049   | STIMEX-TEX 1-20, 1 connector<br>(approx. 23-30 cm ø)        | 2        |
| 107051   | STIMEX-TEX 1-30, 1 connector<br>(approx. 28-40 cm ø)        | 2        |
| 107052   | STIMEX-TEX 1-35, 1 connector<br>(approx. 32-55 cm ø)        | 2        |
| 107053   | STIMEX-TEX 2-20, Elbow 2 connectors<br>(approx. 23-30 cm ø) | 2        |
| 107054   | STIMEX-TEX 2-35, Knee 2 connectors<br>(approx. 32-55 cm ø)  | 2        |

Delivery Content

| Item no | Description                              | Quantity |
|---------|------------------------------------------|----------|
| 104064  | EMP4 ECO +                               | 1        |
| 106351  | Cable type 5.15 (black)                  | 1        |
| 106352  | Cable type 5.16 (grey)                   | 1        |
| 106353  | Cable type 5.17 (red)                    | 1        |
| 106354  | Cable type 5.18 (yellow)                 | 1        |
| 101062  | Charger for EMP4 ECO +                   | 1        |
| 283400  | STIMEX self-adhesive electrodes 50x50 mm | 2        |
|         | Transport case                           | 1        |
| 101723  | Instruction manual                       |          |

## ELECTRODES POSITION REHABILITATION PROGRAMS

P11  
P12  
P13  
P14  
P15  
P16


Hand flexor

Hand extensor

Biceps


Trapezius

Back

Latissimus

Stomach muscles

Gluteal muscles

## ELECTRODES POSITION REHABILITATION, COMBINATION PROGRAMS

P21

P22

P23

P24


P25

Leg extensor

Leg flexor

Foot flexor

Foot extensor

## ELECTRODES POSITION DYNAMIC PROGRAMS

Pain Treatment with  
dynamic stimulation

P1


P3

P9


P20

Leg flexor


Foot flexor (calf)


Leg flexor and foot flexor


Back

Leg extensor (quadriceps)

Back


Back and quadriceps

ELECTRODES POSITION  
PAIN TREATMENT  
PROGRAMS


- P2
- P4
- P5
- P6
- P7
- P8
- P10


Neck

Trapezius


Deltoid


Elbow

Hip

Knee


## ELECTRODES POSITION PAIN TREATMENT PROGRAMS

P2  
P4  
P5  
P6  
P7  
P8  
P10


Lumbosciatica

Ankle

Wrist



## ELECTRODES POSITION INCONTINENCE PROGRAMS

P17

P18

P19

## Program Overview

|                  |                                                          |
|------------------|----------------------------------------------------------|
| Program P1 + U1  | Low Back Pain                                            |
| Program P2 + U2  | Lumbosciatica                                            |
| Program P3 + U3  | Cervicobrachial Syndrome                                 |
| Program 4 + U4   | Chronic Pain                                             |
| Program 5 + U5   | Gate Control + Endorphin Release                         |
| Program 6 + U6   | Gonarthrosis / Coxarthrosis                              |
| Program 7 + U7   | Epicondylitis                                            |
| Program 8 + U8   | Algodystrophy                                            |
| Program 9 + U9   | Anti-Adaptation – Dynamic Stimulation                    |
| Program 10 + U10 | Classic Gate Control 120 Hz                              |
| Program 11 + U11 | Treatment of Contractures                                |
| Program 12 + U12 | Atrophy Upper Extremities                                |
| Program 13 + U13 | Atrophy Lower Extremities                                |
| Program 14 + U14 | Muscle Strengthening Upper Extremities                   |
| Program 15 + U15 | Muscle Strengthening Lower Extremities                   |
| Program 16 + U16 | Venous Reflux                                            |
| Program 17 + U17 | Urge Incontinence                                        |
| Program 18 + U18 | Mixed Incontinence                                       |
| Program 19 + U19 | Stress Incontinence                                      |
| Program 20 + U20 | Relaxation – Dynamic Stimulation                         |
| Program 21 + U21 | Combination: Active wrist joint motion + pain management |
| Program 22 + U22 | Combination: Active elbow joint motion + pain management |
| Program 23 + U23 | Combination: Active knee joint motion + pain management  |
| Program 24 + U24 | Combination: Active ankle joint motion + pain management |
| Program 25 + U25 | Combination: Hemiparesis (shoulder) + pain management    |

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Wuxi Jiajian Medical Instruments Co., Ltd,