## Vigorimeter

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Directions for Use

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Mode d'emploi

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V. 1.1-DE-EN-ES-FR-IT (06.12)

**KLS martin**  
GROUP

**Inhaltsverzeichnis****Table of Contents****Índice****Sommaire****Indice**

<b>Deutsch</b>	<b>Seite 2</b>
<b>English</b>	<b>Page 8</b>
<b>Español</b>	<b>Página 14</b>
<b>Français</b>	<b>Page 20</b>
<b>Italiano</b>	<b>Pagina 26</b>

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**Table of Contents – English**

1	Description .....	9
1.1	Components / Overview of the KLS Martin Vigorimeter .....	9
1.2	Rubber Ball Material Characteristics.....	10
1.3	Adjustment of the Pressure Gauge .....	10
1.4	Information regarding Cleaning and Reprocessing .....	10
1.4.1	Pressure gauge and tube with plug cone .....	10
1.4.2	Rubber balls .....	10
2	Measuring method.....	11
3	Simulation during test.....	11
4	Finger testing .....	11
5	Standard values for healthy people.....	12
6	Range of application .....	13
6.1	Central nerve troubles .....	13
6.2	Troubles mainly caused by diseases of the spinal cord .....	13
6.3	Troubles of the muscular system.....	13

## 1 Description

In the diagnostic field, the testing of the muscular strength plays a most remarkable part.

The Vigorimeter from KLS Martin is a quite new dynamometer the results of which are registered by a pointer on the dial of a manometer. This manometer is connected with rubber bulbs of various sizes by means of a rubber tubing with metal adapter.

The essential advantage of the Vigorimeter from KLS Martin lies in the fact, that it can be easily handled, that the results are absolutely correct and in its wide range of application, that is to say, from the baby to the adolescent and the adult.

The test does not entail any inconvenience and may be, therefore, repeated. Even if the tests are being repeated several times, there will be no troublesome accompanying symptoms.

Owing to the fact that the size of the rubber bulb (small, medium and large) adapts itself to the individual size of the hand the patient is enabled to employ his maximum strength.

The large bulb will be used by adults and adolescents whose metacarpal circumference is more than 19 cm (glove size more than 7 1/2), the medium size bulb by young ladies and women as well as by adolescents whose metacarpal circumference is less than 19 cm, while the small bulb should be used in testing babies and schoolboys and schoolgirls of from 2 to 10 years approximately. Further, this small bulb can be also used for testing the strength of individual fingers, that is to say, of adolescents and adults. In this case, the bulb will be held between the thumb and the fingers or between the metacarpus and the finger.

### 1.1 Components / Overview of the KLS Martin Vigorimeter

#### The Vigorimeter 17-400-01 is composed of:

Pressure gauge	09-400-01-04
Rubber ball large	17-400-95-98
Rubber ball medium size	17-400-90-98
Rubber ball small	17-400-94-98
Cone insert (1 piece each) for the rubber balls	17-400-92-98
Tube with plug cone	17-400-91-98
Spare rubber balls (set of 3)	17-400-05-04

## 1.2 Rubber Ball Material Characteristics

- The above-mentioned item consists of soft PVC (medical-grade Soft-PVC, PVC-P in short),
- the item is latex-free,
- is physiologically harmless as recommended by XLVII-BGA,
- is white, lead-free and cadmium-free,
- does not contain the plasticizer DEHP,
- and is CE-marked.

## 1.3 Adjustment of the Pressure Gauge

Before shipment, the pressure gauge is factory-tested with a very exact and pre-defined test pressure in a so-called 3-point calibration. This means unpressurized, with the mean value and with the highest value (end value). This way, the pressure gauges are adjusted to these values.

“Normal” environmental conditions and intentional use have no influence on the pressure gauges. Mechanical influences (falling, vibrations, etc.) strong temperature fluctuations, humid or also dusty environmental conditions can affect the accuracy of the pressure gauge.

Wrong or inaccurate readings can also be caused by defective tubes or balls.

In order to ensure steady measurement and readings, the pressure gauges must regularly (every 1 – 2 years) be sent in to Gebrüder Martin for checking purposes, cleaning (inside of the pressure gauge) and re-adjustment (as described above).

## 1.4 Information regarding Cleaning and Reprocessing

### 1.4.1 Pressure gauge and tube with plug cone

The pressure gauge (09-400-01-04) and the tube can be cleaned on the outside using a damp but not wet cloth and a common cleaning and disinfection agent.

### 1.4.2 Rubber balls

- Sterilization by boiling: Dismount balls, immerse in boiling water and boil approx. 10 minutes
- Gas sterilization: Sterilization by means of ethylene oxide is possible
- Radiation sterilization: Gamma-ray sterilization is possible (depending on ball color there may be slight color changes)
- Steam sterilization: **not possible**

## **2 Measuring method**

The muscle force is transferred from the hand to the rubber ball, and the air pressure from the rubber ball to the pressure gauge. The ball is tightly fitted on the hose and pressure gauge. Both parts of the cone are tightly fitted together by a quick turn so that air cannot escape.

The patient is given the ball that has already been connected to the pressure gauge in a way that the ball rests in the patient's hollow hand with its largest volume and the hose on the neck of the ball protrudes between thumb and index finger. The fingers that surround the ball should be quickly closed to a fist and by doing so should press the ball together. The finger tips should also be involved when closing the fingers to a fist, that means they should firmly press into the ball.

## **3 Simulation during test**

By using the bulb of medium size, it is possible to employ a greater strength than with the large size bulb. When comparing the results obtained with both these bulbs, you arrive at a difference of 20 – 30 %. However, from the physical point of view, positive rules or relations cannot be established.

In case where by no obvious reason the values found should prove to be extremely low, it will be advisable to repeat the test, using another bulb. If the result should considerably deviate from the difference of 20 – 30 %, the patient should be suspected of simulation, an exception being isolated lesions of some peripheric arm nerves. Physiologically taken, an insignificant decline may happen (less than 10 %) in case measuring is being successively performed without any interval with the same bulb in the same hand. However, after some minutes rest, this decline will be entirely compensated. In case of several repeated tests, the results will rather often be better.

## **4 Finger testing**

The patient should exert a slight pressure between the thumb and the respective finger, while the practitioner will hold the neck of the bulb.

Another possibility:

Place the small bulb under the prominence of the respective finger, the bulb pointing at the wrist (the thumb will not be used). The pressure exerted should be as strong as possible. Test the thumb as follow: Place the small bulb between thumb and metacarpus II (finger outstretched) and perform an energetic adduction. In order to test the force of the opposition, the patient should place the bulb between the prominences of the fingers IV and V, while the thumb exerts an energetic pressure on the bulb.

## 5 Standard values for healthy people

Patient	Ball size	Measurement of hand force	Finger check	Thumb check
				
Men 	Large	0,8 to 1,3 bar	0,05 to 0,4 bar	0,1 to 0,5 bar
Women 	medium	0,7 to 1,25 bar	0,05 to 0,3 bar	0,1 to 0,3 bar
Adolescents 	medium	0,8 to 1,1 bar	0,05 to 0,3 bar	0,1 to 0,2 bar
Children 	Small	0,1 to 0,4 bar	0,05 to 0,2 bar	0,05 to 0,15 bar

## **6 Range of application**

Measuring of the muscular strength is indispensable in diagnostinating the following diseases:

### **6.1 Central nerve troubles**

- a) Traumatic cerebral lesions (Contusio) accompanied by spasms or by latent pareses, postoperative crises after cerebral interventions,
- b) cerebro-organic processes as a consequence of inflammations, encephalitis, encephalomyelitis, tumors, embolisms or mollifications, Little's disease, mobility disturbances as a consequence of a disease of the extrapyramidal system.

### **6.2 Troubles mainly caused by diseases of the spinal cord**

- a) Traumatic cervical and bulbar disturbances, contusio medulæ, paraplegic lesions,
- b) pictures of the general condition after inflammatory diseases, myelitis, acute ascending spinal paralysis (Landry's paralysis) tumors, syringomyelia, troubles of the spinal column in connection with internal diseases, intoxications,
- c) diseases in connection with the material spinal systems, spastic spinal paralysis, progressive spinal muscular atrophy, amyotrophic lateral sclerosis, neural muscular atrophy, myopathy.

### **6.3 Troubles of the muscular system**

- a) Traumatic: muscular laceration, crises after fractures, atrophy by disuse.
- b) muscular diseases, especially progressive muscular dystrophy (ERB), myasthenia.

This enumeration of indications is not complete, however, it shows that thanks to its extraordinary method of application a dynamometer of the type as above is most appropriate to control the patient's rehabilitation, or to help him to recover his muscular strength.

Therefore, this is a diagnostic means which will be welcomed by the majority of all practitioners.

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**Stuckenbrock Medizintechnik GmbH**
**Martin Nederland/Marned B.V.**
1270 AG Huizen · The Netherlands
**Gebrüder Martin GmbH & Co. KG**
**Rudolf Buck GmbH**
**KLS Martin UK Ltd.**