SCALA RIDOTTA "PONTICELLO" 260381

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ESCALIER RÉDUIT "PONTICELLO" 260381

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REDUCED LENGTH WALKING RAMP 260381

Instruction for the assembly

VERKÜRZTE TRAININGSTREPPE: 260381

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CLASSIFICATION AND USE

REDUCED LENGTH WALKING RAMP 260381

Read the instructions and use the product as directed according to its function.

Do not use the product in a way other than its intended purpose.

INTENDED USE

Stairs for walking rehabilitation and muscle reinforcement with possibility to apply an inclined platform. **INTENDED FOR USE IN SMALL AREAS.**

Patients who are able to walk with parallel bars and need to reinforce their muscles as well as loosen their lower limbs. The platform supplied allows to simulate the degree of difficulty of an inclined path.

TYPE OF PATIENT

Patients who have achieved an advanced level of motor recovery - where surveillance by a professional remains indispensable.

PLACE OF USE

Hospitals, specialised clinics, outpatient clinics.

OPERATORS QUALIFIED TO USE THE PRODUCT

All healthcare professionals, physiotherapists, and related industries.

The instructions given to the patient by the reference healthcare professional suffice.

CONTRAINDICATIONS

Follow the safety instructions in user manual.
Not suitable for sports use.

OTHER INDICATIONS

The device does not have essential levels of performance.

**Report all serious incidents that occur due to the device.
Report the incident to the manufacturer and the competent authority**

At the end of its life, the product must be disposed of in accordance with local regulations and taken to specific waste collection facilities.

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Weight: 97 kg

Storage: dry, dark place, -20°+60°

Cleaning: it is generally recommended to clean the device after each use and to carry out the disinfection cycle periodically. To clean the device, wipe its parts with a cloth moistened with warm water. Do not use any agents that could damage the surface of the device through abrasion or overly aggressive treatment.

Operating temperature: 0° / +40°

Expected lifetime: it is estimated that the device can be used for 10 years.

It is mandatory for the user and the patient to use disposable gloves and shoes with clean rubber soles with disposable overshoes

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info@ferrox.it

SRN: IT-MF-000011243

Under its sole responsibility state that the products:

Name	Code	UDI-DI Basic:
"PONTICELLO" SHORT STAIRWAY REDUCED	26038I	805347043SCALA9R
WALKING RAMP	260385	805347043SCALA9R

meet all the applicable essential requirements, of Annex I, implemented by the
(UE) REGULATION 2017/745 OF THE EUROPEAN PARLAMENT AND COUNCIL on Apr 5th 2017.

Classified (Annex VIII, MDR 2017/745):

Class I (Rule 1)

Procedure for the assessment of conformity:

Technical documentation as per annex II and III.

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FERROX S.r.l.